

MIDAS GUI — Configuration & Defaults Guide

A step-by-step guide to setting up your own defaults, paths, materials and algorithms.

Focused how-to companion to the main `gui_documentation.md` (§15). Everything is stored in **one per-user JSON config** that overrides the shipped built-in defaults. You can manage it entirely from the **Preferences dialog**, or edit / share the JSON file directly.

1. What you can configure

Without editing any code:

Domain	Examples
Geometry defaults	wavelength λ , pixel size, Lsd, beam centre
Menus	the clickable pixel-size presets and K-edge foils
Materials	your phases for ring simulation (lattice + space group)
Calibrants	the calibrant list on the Calibrate tab
Default paths	default data / calibration / output files & folders
Algorithms	default calibration pipeline, integration kernel, output format, error model, colormap
Visible tabs	which optional tabs are shown (Data Viewer, Mask, Calibrate, Batch are always on)
Interface scale	whole-application zoom (layout + fonts) for HiDPI / 4K monitors

2. How it works

There is **one per-user config file**. If present, its values override the shipped built-in defaults; if absent, you get the built-ins. There is no separate "group" mechanism — to share settings, export a JSON and hand it to a colleague (§6).

Changes take effect on the next launch (values are read when the window is built). A malformed config file is ignored (built-ins are used) — it will never stop the GUI from starting.

Where the file lives

OS	Path
Linux	<code>~/.config/midas_gui/config.json</code>
macOS	<code>~/Library/Application Support/midas_gui/config.json</code>
Windows	<code>%APPDATA%\midas_gui\config.json</code>

You normally never touch this path by hand — the Preferences dialog reads and writes it. **Settings ▶ Open config folder** jumps straight to it.

3. Method A — the Preferences dialog (recommended)

1. Launch the GUI (`midas-gui` , or `python -m midas_gui`).
2. Open **Settings ▶ Preferences...** (menu bar, top-left).
3. The dialog opens **pre-filled with the full shipped defaults** — you edit from a complete starting point.
Work through the tabs:
 - **Geometry** — wavelength, pixel size, Lsd, beam centre.
 - **Paths** — default files/folders (use the ... button to browse).
 - **Materials and Calibrants** — **Add** a row and fill `name`, `a`, `b`, `c`, `α`, `β`, `γ`, `SG` (`SG` = space-group number), **Remove selected**, or edit any cell. Materials appear in the Data Viewer's material dropdown; calibrants in the Calibrate tab list.
 - **Menus** — the pixel-size presets (label + μm) and K-edge foils (element + keV).
 - **Algorithms** — default calibration pipeline, integration kernel, output format, error model, and colormap/theme.
 - **Tabs** — tick which optional tabs are visible. Data Viewer, Mask, Calibrate and Batch Integrate are always on (locked). This setting **applies immediately** on Save — no restart needed (unlike the others).
 - **Display** — the **interface scale** (whole-app zoom for HiDPI / 4K screens): pick a value or preset (100 % $\approx$ 1080p, 150 % $\approx$ 1440p, 200 % $\approx$ 4K). It applies after a restart; you're offered one on Save when it changes. A quick picker is also on **Settings ▶ Interface scaling....**
4. Handy buttons (top of the dialog):
 - **Save current GUI state** — copies the Data Viewer's *current* λ / pixel / Lsd / beam centre into the Geometry fields (set a system up once, then capture it).
 - **Save config to JSON... / Load config (JSON)...** — export or import a config file (how you share settings — see §6).
 - **Reset to shipped defaults** — discards your config and returns to the built-ins.
5. Click **Save as my defaults**. Your per-user file is written; **restart the GUI** to apply.

Because the tables are pre-filled with the shipped list, "remove" and "modify" just work — whatever is in the tables when you Save becomes your complete list. To get the original shipped list back at any time, use **Reset to shipped defaults**.

4. Method B — edit the JSON file by hand

Prefer the dialog for materials/calibrants (it pre-fills the full list). Hand-editing is convenient for a few scalar overrides.

1. Copy the template to your per-user path (§2), e.g. on Linux:

```
mkdir -p ~/.config/midas_gui
cp documentation/config.example.json ~/.config/midas_gui/config.json
```

2. Keep **only the keys you want to change** — every section and key is optional.
3. Save and (re)launch the GUI.

Format

```
{
  "geometry": {
    "wavelength_A": 0.39, "pixel_um": 75.0, "lsd_um": 121000.0, "bc_y": 10.0, "bc_z":
    10.0,
    "pixel_presets": [["GE", 200.0], ["Varex", 150.0], ["Pilatus", 172.0], ["Eiger",
    75.0]],
    "k_edge_foils": [["W", 69.525], ["Au", 80.725], ["Pb", 88.005]]
  },
  "materials": { "Ni (FCC)":
    {"a":3.5238,"b":3.5238,"c":3.5238,"alpha":90,"beta":90,"gamma":90,"sg":225} },
  "calibrants": { "CeO2":
    {"a":5.4116,"b":5.4116,"c":5.4116,"alpha":90,"beta":90,"gamma":90,"sg":225} },
  "paths": {
    "calibrant_tif": "", "calibrant_h5": "", "nickel_h5": "", "nickel_dir": "",
    "nickel_frame0": "", "calib_file": "", "pdf_iq_file": "", "pdf_calib": ""
  },
  "ui": {
    "calibration_pipeline": "one_shot", "integration_kernel": "subpixel2",
    "output_format": "csv", "azimuthal_method": "poisson", "plot_theme": "hot",
    "visible_tabs": ["Calib. Refinement", "Pump Probe"],
    "ui_scale": 1.0
  }
}
```

Key reference

- **geometry** — numbers in Å / μm / px. `pixel_presets` = [label, size_um] list for the clickable **px** menu; `k_edge_foils` = [element, K-edge keV] list for the clickable **λ** menu.
- **materials / calibrants** — keyed by display name; each is `a`, `b`, `c` (Å), `alpha`, `beta`, `gamma` (deg), `sg` (space-group number). Paths accept `~` and `$ENV_VARS`.
- **paths** — override the default file/folder each tab opens with. Empty `""` (or omit) keeps the built-in.
- **ui** — `calibration_pipeline` ∈ `one_shot` | `first_time` | `four_stage` | `bayesian` | `joint`; `integration_kernel` ∈ `hard` | `subpixel2` | `subpixel4` | `polygon`; `output_format` ∈ `csv` | `xye` | `fxye` | `dat` | `h5` | `2d_csv`; `azimuthal_method` (error model) ∈ `poisson` | `azimuthal` | `hybrid`; `plot_theme` (colormap) ∈ `hot` | `gray` | `viridis` | `inferno` | `plasma` | `turbo`. `visible_tabs` = the list of **optional** tabs to show (any of `Calib. Refinement`, `Corrections`, `PDF Analysis`, `Texture`, `Pump Probe`, `Results & Export`); the four always-on tabs are implicit. The shipped default is `["Calib. Refinement", "Pump Probe"]` (the other four ship hidden). Unlike the

rest, this one applies immediately (no restart). `ui_scale` = whole-interface zoom (0.5–4.0) applied at startup via `QT_SCALE_FACTOR` (~1.5 for 1440p, ~2.0 for 4K); takes effect on restart.

Important — lists replace. If you include `materials`, `calibrants`, `pixel_presets`, or `k_edge_foils`, that section becomes your **complete** list (it replaces the built-in one). Omit a section to keep the built-in list, or use the Preferences dialog, which pre-fills the full shipped list so you always start complete. Valid JSON only — no comments, no trailing commas.

5. Reset to the shipped defaults

- **In the dialog:** Preferences ▶ **Reset to shipped defaults** (then restart).
- **By hand:** delete your per-user `config.json` (§2) and relaunch.

Either returns every setting to the built-in shipped values.

6. Sharing settings with colleagues

There is no shared-file mechanism to configure — you just move a JSON around:

1. Set everything up via the Preferences dialog.
2. **Preferences ▶ Save config to JSON...** and choose a location (email it, drop it on a share, commit it to a project repo — your choice).
3. A colleague opens **Preferences ▶ Load config (JSON)...**, picks the file, reviews it, and clicks **Save as my defaults**.

(Or simply copy the file into their per-user config path from §2.)

7. Common recipes

- **Set a beamline's default energy and detector:** Preferences ▶ Geometry set `wavelength_A` and `pixel_um` (or click the **px** label → your detector) → Save → **Save config to JSON...** and share it.
- **Add a house sample phase:** Preferences ▶ Materials ▶ Add → lattice + SG.
- **Default to polygon integration + XYE output:** Preferences ▶ Algorithms → kernel `polygon`, output `xye`.
- **Point at your data folders:** Preferences ▶ Paths → set the sample/calibration files & folders.

8. Troubleshooting

- **My change didn't apply.** Restart the GUI (settings load at startup).
- **The GUI ignored my file.** It's probably invalid JSON — the app falls back to built-ins. Re-save through the Preferences dialog (it always writes valid JSON).

- **I lost the shipped materials/calibrants.** A hand-edited file with a partial `materials` (or `calibrants`) list replaces the built-ins. Open Preferences (or **Reset to shipped defaults**) to get the full list back, then re-curate.
- **Back to factory defaults.** Preferences ▶ **Reset to shipped defaults** (or delete the per-user `config.json`).